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NOT / NOTES

Usage Instructions 26
Meaning of Symbols 26
Responsibility 26
CE Declaration of Conformity 26
Intended Use 26
Unauthorized Use 27
Safety Warnings 27
Product-Specific Safety Warnings 28
Opening the Box 32
Contents of the Box 32
Technical Information 32
Product Parts Introduction 33
Display Description 33
Indicator Definitions 34
Setup Instructions 35
How to Use 38
Device Cleaning and Maintenance 43
Device Storage 44
Recycling 44
Troubleshooting 44
Shipping 45

Dear Customer,

We hope you enjoy using your product. Thank you for choosing **Fakir Hausgeräte**, a leader in innovative technologies. Please read this user manual carefully before using the product and keep it for future reference.

Usage Instructions

These usage instructions are prepared to ensure that the device is used properly and safely. For your safety and the safety of others, please read the user manual before operating the device.

Meaning of Symbols

In this user manual, you will see the following symbols:

WARNING!

This symbol indicates hazards that could lead to injury or death.

CAUTION!

This symbol indicates hazards that could cause damage to your product or malfunction.

NOTE!

Indicates important and useful information in the user manual.

Responsibility

Fakir Elektrikli Ev Aletleri San. ve Tic. A.Ş. accepts no responsibility for any damage or repairs made outside Fakir Authorized Service Centers or for improper use of the device resulting from non-compliance with the instructions in this user manual.

CE Declaration of Conformity

This device complies fully with the European Union's 2014/30/EU Electromagnetic Compatibility and 2014/35/EU Low Voltage Directives. The product bears the CE mark on the device type label.

Intended Use

The product model is designed for home use. It is not suitable for workplace or industrial use. Any malfunction resulting from such use will not be covered by our authorized service.

Unauthorized Use

In the following cases of unauthorized use, the device may be damaged or cause injuries:

- **Restricted Capabilities and Supervision:** This device should not be used by persons (including children) with limited physical, sensory, or mental capabilities or lack of experience and/or knowledge, unless they are supervised by a person responsible for their safety or given instructions regarding the safe use of the device. Children should be supervised to ensure they do not play with the device.
- **Use with Limited Abilities:** The device can be used by children aged 18 years and above and persons with limited physical, sensory, or mental capabilities or lack of experience and knowledge only if they are supervised, given instructions on safe use, and understand the potential hazards.
- **Packaging Materials:** Keep materials like bags included in the package away from children. There is a risk of suffocation!
- **Cleaning and Maintenance:** Children under 18 years should not clean or maintain the device unless they are supervised.
- **Device and Cable Safety:** Keep the device and its cable out of the reach of children under 18 years old. Children should be supervised to ensure they do not play with the device. Electrical devices are not toys for children. Therefore, use and store the device out of the reach of children. Children may not be aware of the dangers associated with the use of an electrical device. Ensure the power cord does not hang down from the device.

Safety Warnings

When using an electrical device, always follow the precautions below to prevent fire, electric shock, and personal injury:

- **Read the User Manual:** Carefully read the information in the user manual and keep it for future reference. If you give the device to a third party, be sure to include the manual.
- **First Use:** When you remove the product from its packaging, make sure all accessories and attachments are complete. If any parts are missing, contact your nearest Fakir Authorized Service Center. Keep the protective bags of the accessories and attachments away from children and infants

to prevent the risk of suffocation.

- **Check the Condition of the Device:** If you suspect that there is any issue with the device, its casing, cable, or plug, do not operate the device and contact the nearest Fakir Authorized Service Center.
- **Cable and Plug Safety:** Do not use the device if the power cord or plug is faulty or damaged. Do not use the device with an extension cord and ensure that the power cord does not come into contact with heat sources or sharp edges. When unplugging, always pull the plug itself, not the cord. Do not pull on the cord or lift the device by the cord.
- **Child Safety:** Never leave the machine unattended while it is plugged in. Do not allow children to play with the device unsupervised. Store the device's power cord out of the reach of children. Electrical devices are not toys for children.
- **Use of the Device:** Do not use the device outdoors. Your machine may be adversely affected by dirt or rain. Do not immerse the machine in water and protect it from rain or moisture. The device operates with electrical power, which poses an electric shock hazard. Therefore, do not touch the plug with wet hands. Do not use the device to clean flammable or explosive materials (e.g., gasoline or alcohol), burning or smoldering objects (e.g., cigarettes, matches, or hot ashes).
- **Maintenance and Cleaning:** Before cleaning and maintenance, unplug the device and allow it to cool down.
- **Usage Restrictions:** Use the electrical device only as intended. Using the device contrary to the instructions may shorten its lifespan and pose safety risks.

Product-Specific Safety Warnings

- Please read and follow these instructions to optimize device performance and prevent or minimize the risk of injury, skin reactions, and side effects.
- This device is not a toy; keep it out of reach of children.
- Inspect the device before use. DO NOT USE if any damage is detected.
- Keep the device and the power adapter dry.

- Ensure that the operating voltage of the device is compatible with the rated voltage of your power supply.
- Cleanse your skin before use; ensure the skin is completely dry and free from hair, lotions, or gels.

Do Not Use Under the Following Conditions

The following individuals should not use this device:

- Persons under the age of 18 or over the age of 65.
- If your skin tone is marked with an “X” in the “Suitable Skin Tones” table at the beginning of the user manual.
- Individuals who are pregnant, breastfeeding, or menstruating.
- Individuals with active implants such as pacemakers, neurostimulators, or insulin pumps.
- Individuals with sensitive skin prone to rashes or allergic reactions.
- Persons with Skin Type V and VI (very dark skin that rarely sunscreens). In these cases, there is a high risk of hyperpigmentation, hypopigmentation, severe redness, or burns.
- Individuals with a history of skin disease or those currently undergoing dermatological treatment.
- Persons with pre-cancerous lesions or multiple atypical moles (nevi) in the treatment area.
- Persons with infections, burns, eczema, folliculitis, abrasions, open wounds, herpes simplex, cuts, or hematomas in the application area.
- Individuals with a history of collagen disorders, keloid scar formation, or impaired wound healing.
- Persons with vascular diseases such as varicose veins.
- Individuals with diabetes, lupus erythematosus, porphyria, congestive heart failure, or bleeding disorders.
- Individuals with a history of immunosuppressive diseases (including HIV infection and AIDS).
- Persons who are currently receiving or have received radiotherapy or chemotherapy within the last few months.

DO NOT USE the device on the following areas:

- Areas with hyperpigmentation, hypopigmentation, intense redness, or burns.
- Areas where surgical procedures have been performed.
- Over eyebrows or moles that have sustained injury.
- Nipples, areola, labia minora, vagina, anus, and inside the nostrils or ears.
- The scrotum and facial area.
- Over silicone implants, pacemakers, subcutaneous injection ports (e.g., insulin pumps), or body piercings. Avoid contact with any artificial surfaces.
- Over moles, freckles, large veins, dark pigmented areas, scars, and skin abnormalities (without medical supervision).
- Areas with tattoos or permanent makeup; this may cause skin burns and pigment changes.
- Areas where long-lasting deodorants are used; this may lead to skin irritation.
- On sunburned, recently tanned skin, or skin treated in tanning beds (solariums).

DO NOT USE the device if you are taking the following medications or undergoing these treatments:

- If your skin is being treated with Alpha-Hydroxy Acids (AHA), Beta-Hydroxy Acids (BHA), topical isotretinoin, or azelaic acid.
- If you have used oral or topical Accutane or Roaccutane within the last 6 months.
- These treatments make the skin more susceptible to tears, trauma, and irritation.
- If you are taking any medication or substance that causes photosensitivity; please read the medication's package leaflet carefully.
- Do not use the device if there is a risk of light allergy or phototoxic reactions, or if you are advised to avoid sun exposure during medication use.
- If you are using blood-thinning medications (including high-dose aspirin) and cannot discontinue them at least one week before the procedure.

Please Pay Attention to the Following:

- This information is for general guidance and is not exhaustive. If you experience any discomfort during use, discontinue use immediately and consult a healthcare professional.
- This product is not a medical device and is not intended to diagnose, treat, or prevent any disease.
- The performance of the product may be limited on white or grey hair.
- Caution is advised before use on dark skin tones. A test should be performed at a low setting before the first application.
- Avoid direct sun exposure before and during the course of using the device.
- Before the first use, test the device on a small area of your arm at a low power/sensitive mode to observe your skin's reaction
- Ensure that the light exit window is in full contact with the skin during application.
- It is recommended not to use the product in very low-light or dark environments.
- To ensure the longevity of the device, do not operate it under direct sunlight.
- After use, position the device with the ventilation outlet facing up and the application window facing down to help dissipate heat and cool the device.
- Do not subject the device to heavy impacts and do not drop it.
- Turn off the device before unplugging it from the electrical outlet; otherwise, the device may be damaged.

Carefully read the safety rules and usage instructions in this user manual. **Keep this user manual for future reference.**

Opening the Box

- Before first use, remove the device and accessories from the box and check for damage or defects. If the device is damaged or defective, do not use it and contact **Fakir Authorized Service**.
- Repairs should only be performed by Fakir Authorized Service Centers.
- Repairs by unauthorized persons may be harmful to the user.

WARNING!

Do not use a faulty or defective product.

Contents of the Box

• Auria IPL Hair Removal Device	1 pc.
• User Manual and Warranty Certificate	1 pc.
• Storage Pouch	1 pc.

Technical Information


Light Exit Window	: 3.0 cm²
Adapter Input	: 100–240 V AC, 50/60 Hz, 1,5 A Main Unit Input: 24 V DC, 2,7 A 64.8W
Wavelength Range	: 560 ~ 1200 nm
Total Number of Pulses	: 999.999
Maximum Energy Density (Fluence):	
Single Pulse Mode	: 6,3 J/cm²
Double Pulse Mode	:6,4 J/cm²

Pulse Duration:	
Pulse Duration	: < 12 ms
Single Pulse Interval	: 250 ms
Output Energy:	
Single Pulse Mode	: 0,257 ~ 19 J
Double Pulse Mode	: 0,509 ~ 19,2 J

NOTE!


The values indicated on the device or in other printed documents provided with the product are obtained under laboratory conditions according to relevant standards. These values may change depending on the product's use and environmental conditions.

Product Parts Introduction

1. Air Outlet
2. Light Exit Window (Sapphire Cooling Surface)
3. Skin Tone Detection Area
4. Flash (Pulse) Button
5. Air inlet
6. LED Screen
7. Double Pulse On/Off Button
8. Power On/Off and Intensity Control Button
9. Power Input Port
10. Flash (Pulse) Ready Indicator
11. Power Cable
12. Power Adapter

Display Description

1. Intensity Level Indicator

- Indicates the adjustable energy intensity level (Levels 1–5).
- A total of five energy bars is displayed.
- The number of illuminated bars while the device is operating represents the selected intensity level.
- For example, if only one bar is illuminated, the intensity level is set to Level 1.

1.1 SR Indicator (SR)

The SR indicator illuminates when the intensity level is set between Levels 1 and 2. This mode operates at a low energy level.

1.2 HR Indicator (HR)

The HR indicator illuminates when the intensity level is set between Levels 3 and 5. This mode operates at a high energy level.

2. Fan Indicator

- Ensures cooling of the internal components of the main unit while the device is operating.
- When the device is powered on, the fan icon illuminates and the fan starts operating.
- After the device is powered off, the fan continues to run for approximately 15 seconds.

3. Double Pulse

When the indicator is illuminated:

The double pulse mode is active. In this mode, two light pulses are emitted with each trigger.

When the indicator is off:

The single pulse mode is active. In this mode, one light pulse is emitted with each trigger.

4. Number of Pulses Used

- The counter increases by 1 with each light pulse emitted.
- The counter automatically resets when the device is powered off.

Indicator Definitions

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When the green light is steady:

1. In manual flash mode, the attachment is not in contact with the skin.
2. In skin tone detection mode, the attachment is not in contact with the skin.


When the green light is blinking:

1. In manual flash mode, the attachment is in contact with the skin and ready to emit a pulse; however, the 'Flash' button has not yet been pressed.
2. Automatic flash mode is activated.

Setup Instructions


Preparation

- Connect the power adapter to the main unit as shown in the illustration below.
- Plug the adapter into a suitable electrical outlet.
- The device enters standby mode, and “OFF” is displayed on the screen.

Powering On the Device

- Press and hold the Power On/Off - Intensity Control Button for approximately 2 seconds.
- The display screen will illuminate, and an audible beep will sound.
- The device starts operating at the lowest intensity level (Level 1) with double flash mode activated.
- The ventilation fan and the cooling system will activate simultaneously upon startup.

Powering Off the Device

- While the device is on, press and hold the Power On/Off - Intensity Control Button for approximately 2 seconds.
- The screen will turn off, and the "OFF" indicator will be displayed.
- An audible beep will sound, and the device will be turned off.
- The fan continues to operate for approximately 15 seconds after the device is turned off.


Controlling Intensity Level

- While the device is powered on, press the Power On/Off - Intensity Control Button to set the desired intensity level.
- The device features 5 distinct intensity levels (Levels 1–5).
- Each press increases the level by one, cycling back to Level 1 after Level 5.
- First-time users are recommended to start at Level 1.
- The intensity level can be gradually increased according to skin tolerance.


Level 1	Level 2	Level 3	Level 4	Level 5
				

Manual Pulse Mode

- When the device is powered on, Manual Pulse Mode and Double Pulse Mode are automatically activated.
- Place the device probe (treatment head/window) in full contact with the skin area to be treated.

When the Flash Ready Indicator begins to flash, press the "Flash" Button to trigger the light pulse.

Press the “ 

” button to trigger a light pulse.

NOTICE!

Manual Pulse Mode

- If the attachment is not in full contact with the skin, the Flash Ready Indicator illuminates; however, pressing the "Flash" button will not emit a pulse.
- If the attachment is in full contact with the skin but the "Flash" button is not pressed, the Flash Ready Indicator blinks.
- If the light exit window remains in contact with the skin for 8 seconds, the Flash Ready Indicator blinks and the device stop emitting pulses.
- Manual Pulse Mode is recommended for hair removal on small areas.

Automatic Pulse Mode

- While the device is in Manual Pulse Mode, press and hold the "Flash" button for approximately 2 seconds.
- Two short beep sounds will be heard.
- The Flash Ready Indicator starts blinking, and the device switches to Automatic Pulse Mode.
- When the attachment is pressed against the treatment area, the device automatically emits pulses without pressing the "Flash" button.

NOTICE!

Automatic Pulse Mode

- If the attachment is not in full contact with the skin, the Flash Ready Indicator blinks and the device does not emit pulses.
- When the light exit window is in contact with the skin, pulses are continuously emitted for up to 8 seconds; after this period, the device stops emitting pulses.
- Automatic Pulse Mode is recommended for hair removal on large areas.
- Ensure that the entire treatment area is evenly covered during use; do not overlap pulses.

Skin Tone Detection Mode

- Move the device away from the skin and press and hold the "  " button for approximately 2 seconds.
- The level indicator lights will cycle, and the device switches to Skin Tone Detection Mode.
- When the indicator light corresponding to the detected level remains steady, the device will exit the detection mode.
- Prior to starting the treatment, identify your skin tone using the provided Skin Tone Chart.
- The device can be used safely if your skin tone falls within the compatible range specified in the chart.
- When the attachment is pressed against the application area, the device will automatically emit flashes without requiring a manual press of the "Flash" button.

Mode Selection – Manual vs. Automatic Pulse Mode

- While the device is powered on, move the device away from the skin and press and hold the "Flash" button for approximately 2 seconds.
- Two short beep sounds are heard, and the Flash Ready Indicator blinks; the device switches to Automatic Pulse Mode.
- While the device is powered on, move the device away from the skin and press and hold the "Flash" button again for approximately 2 seconds.
- A short beep is heard, and the Flash Ready Indicator remains steady; the device is now in Manual Pulse Mode.

Single / Double Pulse Flash Mode Selection

Press the" " button briefly to switch the flash mode.

If the Double Pulse Indicator is ON:

The device operates in Double Pulse Mode, emitting two pulses per activation.

If the Double Pulse Indicator is OFF:

The device operates in Single Pulse Mode, emitting one pulse per activation.

NOTICE!

When the device is first powered on, the double pulse mode is switched on by default. To disable Double Pulse Mode, press the " " button once.

How to Use

Skin Test

For first-time users, it is recommended to perform a skin test on the arm to observe any reactions before using the device.

a. Turn on the device and select intensity Level 1. Place the attachment in full contact with the skin and apply one pulse.
b. Wait 24 hours and check the tested area for any reactions such as rash, pain, or discomfort.
c. If the skin in the tested area appears normal (no pain or discomfort), you may continue using the device. If any abnormal reaction occurs, discontinue use immediately.
d. Increase the intensity by one level (Do not apply multiple pulses to the same area).

Repeat steps a-c for each intensity level.

NOTICE!

- In Skin Tone Detection Mode, if the attachment is not in full contact with the skin, the device will not emit pulses.
- Skin Tone Detection Mode is recommended for hair removal on large areas.
- In this mode, the device's skin sensor automatically detects skin tone for your safety.
- If your skin tone is not suitable for treatment, the device will not emit pulses.
- If your skin tone falls within the suitable range, the device will emit pulses.
- In Skin Tone Detection Mode, the IPL system re-detects the skin tone each time the light exit window is moved.

Example:

- If the device does not cause any discomfort (or pain) at Level 1, you may proceed to Level 2 following steps a–c outlined above.
- After completing the skin test, use the intensity levels that did not cause any adverse reaction on your skin.
- After each test, wait 24 hours to verify that your skin does not react negatively to the treatment.
- If no visible reaction occurs, you may continue using the device.

Usage Recommendations

Before using the device, if hair has been removed by plucking, waxing, or other permanent methods, wait until hair regrowth begins before using the device.

Treatment Areas

- Small areas, such as arms and underarms.
- Facial hair below the jawline
- Large areas on the body such as legs

Treatment Frequency

1. Initial Phase (First Treatment Cycle)

- Weeks 1–4: Apply once per week.
- Weeks 5–8: Apply once every 2 weeks.
- Weeks 9–12: Apply once every 3 weeks.
- This phase is necessary to target hairs in the active growth (anagen) phase.
- Visible hair reduction may begin after several treatments, depending on skin type and hair structure.

NOTICE!

- Do not perform treatments more frequently than the recommended intervals. Application on the same area at intervals shorter than 7 days is not recommended.

2. Maintenance Phase (Follow-up Treatments)

- After the desired hair reduction has been achieved, it is recommended to perform one application every 4–8 weeks to maintain the results
- The treatment frequency may be adjusted according to hair regrowth and individual needs.

NOTICE!

- Do not treat the same area more than once per day.
- Increasing the frequency of use does not improve effectiveness. Best results are achieved with regular use according to the recommended schedule.
- Results may vary from person to person due to hormonal changes and individual differences.

Skin Tone / Hair Color Reference Chart

Before starting treatment, consult the table below to verify the suitability of your skin tone and hair color for device use.

 : The device is suitable for you.

 : The device is not suitable for you.

Skin Tone Information

- This product is suitable for Skin Types I – IV.
- Please identify your skin tone and consult the table below to check if the device is suitable for you.

Skin Tone Table

Skin Tone/Color	Appearance	Suitability Status
Type I	Fair	Suitable
Type II	Light	Suitable
Type 3	Beige	Suitable
Type IV	Medium brown	Suitable
Type V	Dark brown	Not Suitable
Type VI	Dark	Not Suitable

About Dark Skin Tones

- Darker skin tones absorb light more easily, which may lead to skin irritation.
- Therefore, users with darker skin tones are advised to use a lower flash intensity.

For Darker Skin Tones


- Start at the lowest intensity level.
- Adjust the intensity based on your individual pain tolerance and skin sensitivity.

Hair Color Compatibility

- The IPL system is effective on natural black or brown hair.
- It may have reduced or no effect on white, gray, blonde, or red hair.

Steps

Suitable Hair Colors

Black / Brown / Dark Blonde: Suitable

Light Blonde / Red / White / Grey: Not suitable.

Prepare for Treatment

- Before using the device, trim all hair on the skin surface using a razor or an electric shaver.
- Do not use tweezers, waxing, or chemical depilatory products (Figure 9).
- Thoroughly cleanse and dry the skin.
- Ensure the skin is free from any residues of body sprays, lotions, perfumes, sunscreens, or other cosmetic products.

How to Operate the Device

- Connect the power adapter to the main unit and plug it into an electrical outlet.
- Press and hold the Power On-Off/Intensity Control Button " ⏻ " to turn on the device.
- Select the appropriate flash intensity level based on your skin tone and hair color.
- First-time users are recommended to start at the lowest intensity level and increase the intensity gradually.

NOTICE!


Do not look directly into the light source.

Application

- Press the treatment head firmly against the skin, ensuring full contact between the window and the skin surface.
- Wait for a few seconds; once the Ready-to-flash light starts flashing, press the 'Flash' button to initiate the treatment.

NOTICE!

- **Full Coverage (Correct):** Place the flashes side-by-side like puzzle pieces. Start the next flash exactly where the previous one ended.
- **Overlapping (Incorrect):** Do not overlap the flashes. Overlapping may cause excessive heat exposure on the skin.
- **Missed Areas (Incorrect):** Ensure there are no gaps between flashes. Any areas that are skipped will remain untreated.

After Use

- Once the application is complete, power off the device.
- Unplug the power adapter from the electrical outlet.
- Clean the device and store it in a safe place.

NOTICE!

- This device features a built-in cooling care function.
- It activates automatically when the device is turned on and deactivates when the device is powered off.
- The semiconductor cooling plate is in direct contact with the Sapphire Cooling Plate.
- When the light exit window makes direct contact with the skin, a cooling sensation is felt because the surface temperature of the sapphire plate is lower than the skin temperature.
- This effect is based on the principle of physical heat transfer.
- The device is equipped with an internal fan to cool the internal components of the main unit during operation.
- The fan activates automatically upon startup.
- Once the device is turned on, the fan operates with a distinct airflow sound and continues to run for 15 seconds after the device is powered off to ensure proper cooling.
- After unplugging the power adapter, wait 20 minutes before cleaning the device in accordance with Section 6 of this manual.
- Following the cleaning process, place the device and adapter back into their original packaging and store them according to the specified storage conditions.

CAUTION!

- After using the device, do not apply any scented skincare or cosmetic products to the treated area.
- If the device is used on the underarm area, do not apply deodorant immediately.
- Wait until any skin redness has completely subsided.
- After each application, ensure that no redness persists on your skin for a period of 24 hours.
- Protect the treated area from direct sunlight for 2 weeks following the application.

Device Cleaning and Maintenance

Cleaning
After each use, clean the device according to the following instructions.

Purpose

- This device is designed for the repetitive removal of unwanted hair.
- During use, the device surface may come into contact with organic contaminants and microorganisms carried by hands.
- Therefore, cleaning after each use is mandatory.

Cleaning Method:

- Clean the device immediately after each use.
- Wipe with a slightly damp, soft, lint-free cloth.

Wiping Durations:

- General surfaces: minimum 15 seconds.
- Light Exit Window, maximum 3 seconds.

NOTICE!

a. Ensure that no liquid enters the device during cleaning.
b. Do not use abrasive cleaning agents, thinners, or benzene.
c. Do not immerse the device in water or any other cleaning liquid.
d. Wait at least 10 minutes after cleaning for the device to dry completely.

Maintenance

- The Light Exit Window is a highly delicate component. It is extremely important to protect it from contamination and damage.
- Do not touch, press on, or allow any tools to come into contact with the Light Exit Window.
- This device is not waterproof.
- Do not immerse it in water or wash it under running water.

- Do not use steel wool, abrasive cleaners, chemicals, or corrosive liquids to clean the device.

Device Storage

- If you do not plan to use the appliance for a long while, store it in a cool and dry place.
- Ensure that the appliance is unplugged and its parts are completely dry.
- Ensure that the appliance is kept out of reach of children.

Recycling

When your product reaches the end of its life, follow these steps for proper recycling:

- Cut the device's cable to make it unusable.
- Dispose of the device in designated waste bins for such devices according to the applicable laws in your country.
- Electrical waste should not be disposed of with household waste.
- Throw the device into bins designated for electrical waste.
- The device's packaging is made from recyclable materials.
- Dispose of packaging materials in appropriate recycling bins.
- The device complies with the WEEE (Waste Electrical and Electronic Equipment) Directive.

NOTE!

- Check if you can fix the problem or defect yourself before calling customer service. For this, please refer to the "**Troubleshooting**" section. If you cannot resolve the issue or defect, call customer service.
- Before calling, please have the **model name** and product **serial number** ready. You can find this information on the **product's type label**.

 Help Hotline: (+90 212) 444 0 241

Troubleshooting

Your device has undergone quality control to avoid any defects. However, if the device does not work for any reason, please try the solutions listed in the following table before contacting Fakir Authorized Service.

Problem	Probable Cause	Solution
Screen light does not turn on.	Power cable is not connected correctly.	Check the power adapter connection.
	Power on-off/intensity control button was only pressed briefly.	Press and hold the power on-off/intensity control button for approximately 2 seconds.
Device does not flash.	Power cable is not connected correctly.	Check the power connection.
	The attachment is not in close enough contact with the skin.	
Device enters Standby Mode during use.	No light pulse was emitted for 3 minutes.	Press and hold the power on-off/intensity control button for 2 seconds to restart the device.
Unpleasant odor during use.	The Light Exit Window filter is dirty.	Clean the window with 70% isopropyl alcohol.
	Hair may not have been shaved properly; the flash may be burning remaining hair.	Ensure hair is trimmed with a razor or shaver before treatment.
	Flash intensity is set too high.	Decrease the intensity level.
Skin redness after use.	Mild redness is a normal and harmless reaction to the hair removal process.	No action is required; it will subside naturally.

Shipping

To prevent damage to the device, follow these steps to transport your device safely:

- Transport the device in its original packaging.
- If the original packaging is not available, use a suitable and soft packaging that protects the device.

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Scan the QR code for the user manual in different languages.

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